



# AOX-400



Automatic Organic Halogen Analyzer



Fully automatic solution  
for Analysis of adsorbable Organic Halogens





# Accurate Objective eXtraordinary

Nittoseiko Analytech, former Mitsubishi Chemical Analytech, provides a high-end instrument of the newest generation for automated measurement of Adsorbable Organically bound Halogens: AOX-400.

The screening of AOX sum parameter is quite important to assess the quality of water, sludge or soil for contamination by organic halides (chlorine, bromine, iodine) which are adsorbed to activated carbon. Possible toxic effect from AOX on humans, animals, plants and microorganisms can be prevented by monitoring with AOX-400.

This instrument opens up an easy way to perform environmental analysis by combustion in an open top furnace which guarantees easy maintenance. The usage of air as combustion gas is a unique feature which allows low running costs.

## Measurement principle


AOX adsorbed to activated carbon is pyrolyzed in a quartz combustion tube. Combustion gases containing hydrogen halides are forced through a dehydrating tube (sulfuric acid) by carrier gas and subsequently adsorbed within the electrolyte of a coulometric titration cell.

Within the titration cell the halogen ions are quantified by microcoulometry according to Faraday's law. When hydrogen halide is introduced into the titration cell, the potential changes.



The amount of halogen is calculated from the quantity of electricity required for the generation of silver ions until the end potential is reached.

The gas is moved by application of a suction pump. This principle is used to destroy acetic acid vapor by returning the gases through a separate combustion tube within the furnace after passage of the titration cell.



# High flexibility

Advanced Autosamplers provide flexible combination of different measurement modes easily. AOX-400 offers convertibility between three different injection modes for:

## Shaking method



## Main unit



## Column method



## EOX



By the exchangeable head of sample introduction system SI-400 the measurement of AOX columns, AOX frits and EOX liquid organics in one unit becomes reality. Switching between column injector, frit injector or liquid injector can be performed within a few minutes. Also, the sample tray can be exchanged very quickly. Automatic and direct sample introduction by SI-400 is performed without source of contaminations which opens up low running costs and long usage times.

In addition to :

- **30 Frits**

SI-400 offers the opportunity to measure up to:


- **60 Columns**
- **40 2ml Vials**

which ensures high sample throughput.




# *Favourable design*

AOX-400 is, due to its space saving design, the ideal solution for every lab. Regardless you want to conduct AOX analysis in huge contract laboratory or in a small environmental authority. The small footprint of 46 x 42 cm makes AOX-400 to a compact and easy to handle instrument for environmental analysis.



## Combustion concept

The type of construction makes AOX-400 feasible to move quickly. In addition, maintenance by checking mechanical or electrical components is quite simple. Furthermore, the open top furnace concept provides easy maintenance of the vertical arranged combustion tube to assure long lifetimes. The combustion concept itself is associate with high safety standards: Combustion gases are pulled through the system by a vacuum pump which makes gas leakage almost impossible. Also, the appearance of acid suck-back into furnace can be excluded. With the deodorization tube which is also included in the furnace the vapor of the acetic acid is also combusted. Therefore no unhealthy vapor can enter to the lab air during the operation. Another unique feature of our combustion concept is that no gases except from air are required. Thus, independent operation under cost efficient and sustainable conditions is possible by taking into consideration that 50-75% of initial investment costs additionally spent on gases. The combustion concept is the excellent result of decades of development and is tailor made for complete combustion of all kinds of samples for AOX analysis.



# High precision and reliability





Analytical instruments from Nittoseiko Analytech convince with their precise and accurate measurement results. During instrument development we are focussing on stable and also highly reliable measurement values. One of our strength is the operation under extremely low background levels. Due to vertical combustion, sample boats are not necessary and samples are protected against ambient air by a cover. All in all, AOX-400 presents best performance within a huge measurement range from trace analysis to percentage level.

With AOX-400 argentometric process was increased significantly by titration cell cooling which results also in longer electrolyte lifetime and improved running costs (reagents and consumables). This also leads to less down time.

## Performance

Initial validation data according to ISO official method (9562) show that AOX-400 performs with high linearity and good recovery within the specification of standard methods.

Measurement	Mode	Slope	Correlation coefficient
AOX	Column	0.98	0.99998
EOX	Liquid	0.99	0.99992



※ Reference data using standard sample



# Fastest Analyzer in the market

What makes AOX-400 unique in the market for environmental analysis is the combination of open top furnace concept and automatic gas control. With this particularly feature combination, AOX-400 manages the tough challenge to deal on the one hand with soluble AOX samples and on the other hand with liquid EOX samples effortlessly. Furthermore with AOX-400 it is possible to analyze environmental samples within a very short period of time under high system stability. Due to many optimization regarding automation and titration control AOX-400 makes it possible to measure samples up to 30% faster than our previous AOX-200. This makes our new instrument to an highly efficient analyzer for AOX and EOX.

Measurement	Mode	Detection time [min]
<b>AOX</b>	Frit	8
	Column	8
<b>EOX</b>	Liquid	7



※ Reference data using standard sample

## Official methods


Measurement can be conducted with the following official methods:

- **Column method:** Adsorbable organic halogens (AOX) in water can be measured by DIN EN 1485 or ISO 9562. Total organic halides (TOX) in water can be measured by EPA 9020A.  
Option: Column Adsorption Module Model TXA-04.
- **Batch method:** Adsorbable organic halogens (AOX) in sludge and sediment can be measured by DIN 38414-S18. Adsorbable organic halogens (AOX) in water can be measured by ISO 9562.  
Option: Batch Adsorption Unit Model SA-200.
- **EOX method:** Total extractable organic halides (EOX) in solids can be measured by DIN 38414-S17 and EPA 9023.

# High precision and reliability

The exclusive new software of AOX-400 ensures user-friendly operation, handling and control of measurements.

For startup and shutdown an easy-to-use navigation bar is implemented which takes the user through the measurement. By intuitive shortcut buttons the operation of AOX-400 is self-explanatory so that the user can start the analysis without long training periods. For customized analysis, several different parameters can be adjusted within the software. Control of the measurement can be performed by live view.



The performance of AOX-400 instrument and the titration cell can be checked by automatic quality control functions. For integration of measurements into daily business at laboratory standby functions are available. After the measurement, the results can be reported in different settings.

## AOX-400 at a Glance

- Fastest AOX analyzer in the market
- High sample throughput
- Flexible 3-injection modes are available

**Column method**

**Shaking method**

**EOX**

- Easy maintenance by open top furnace concept
- Low running costs due to unique usage of air as combustion gas
- Unique suction system
- Stable and highly precise measurement results
- Space-saving design
- User-friendly and intuitive software solution



# Specifications

Model	Automatic Organic Halogen Analyzer Model AOX-400
Measuring Object	Ground water, surface water, drinking water, waste water, soil, sludge
Sample forms	Solid (AOX: activated carbon adsorbed liquid sample), Liquid (EOX: extraction solution)
Analysis method	Oxidative pyrolysis / Coulometry
Furnace temperature	300~1000 °C (Regular use: 950 °C)
Sample introduction	Introduction to the open top pyrolysis tube driven by Newton's law AOX: Extrusion injection of activated carbon adsorbed sample (Column method) Drop frit with activated carbon adsorbed sample (Shaking method) EOX: Direct injection by syringe
Measurement range	Total organic halogen: 0.1 to 50 µg
Sample volume	AOX sample: 50 mg or less (1-3x), EOX sample: 200 µL or less
Measurement time	within 10 minutes / measurement (At 2 µg sample measurement)
Number of samples	Frit : 30 samples Column : 60 samples EOX : 40 samples (Size: 2 mL vials with septum)
Gas	Using AIR-200, Air gas is not required. AOX: Air (800 mL / min) EOX: Air (100 mL / min) + Oxygen (300 mL / min)
Operation Condition	15 to 35 °C, 75 %RH or less (No condensation)
Power supply	AC 100/115/230/240 V, 50/60 Hz, 960 VA
Dimensions Weight	AOX-400 main unit : Approx. 460(W) x 420(D) x 560(H) mm, Approx. 37 kg SI-400 : Approx. 200(W) x 305(D) x 410(H) mm, Approx. 7 kg AIR-200 : Approx. 100(W) x 400(D) x 220(H) mm, Approx. 5 kg
Device control and data processing	PC, Printer: compatible with Windows®

# Sampler

Shaking method	Frits set
Column method	SI-400 + Column set
EOX	SI-400 + EOX set

## Nittoseiko Analytech Co., Ltd.

7-10-1 Chuo-rinkan, Yamato, Kanagawa 242-0007, JAPAN  
Tel: +81(0)46-278-0056  
URL: <https://www.n-analytech.co.jp/global/>



Your agent for Switzerland:

**TecDev**  
Le Villaret  
Rte de la Robellaz 6  
CH-1432 Belmont s/Yverdon  
Switzerland

**Phone:** +41 (0) 24 435 22 02  
**Fax:** +41 (0) 24 435 22 03  
**Mobile:** +41 (0) 79 435 96 26  
**Email:** [tecdev@swissonline.ch](mailto:tecdev@swissonline.ch)  
**Visit us:** [www.tecdevexpress.ch](http://www.tecdevexpress.ch)