

Auto Steam Sterilizer

Model: VS-1321 Series

Features

- Microprocessor Control
- Full Automatic Procedure
- Programmability
- Vent Type Auto-Selection
- Appointment Sterilization
- Self Diagnostic
- Water Level Monitor
- Ideal Construction
- Safety Device



Auto-steam Sterilizer

Model	VS-1321-60	VS-1321-80	VS-1321-100
Dimension (In)	Dia.350×630mm	Dia.400×630mm	Dia.450×630mm
(OUT)	670×520×910(H)mm	780×630×940(H)mm	780×630×940(H)mm
Capacity	60 ℓ	80 ℓ	100 ℓ
Heater	2.4Kw	3.6Kw	4Kw
Weight	95Kg	132Kg	144Kg
Power	220V, 11.A, 60Hz	220V, 16.3A, 60Hz	220V, 28.2A, 60Hz

Temp. Range	Ambient to 121 °C
Temp. Accuracy	±1 °C
Operating Pressure	1.0 to 1.2kg/cm ²
Temp. Controller	Microprocessor controller
Display	Back light LCD 10 character display
Program	10 Programs
Timer	0-60min. & Pre-appointment function
Appointment	Pre-appointment function
Pressure Gage	Digital Type, 0-5kg/cm ²
Air Vent	Automatic selectable
Water level monitor	Two point sensor, Auto refill function
Safety Device	Over Heat/ Over pressure protector, water level monitor
Monitor unit	Audible & Visible Device
Basket	Standard Type, 2EA



Autoclave

Autoclave

Model: VS-1221 Series

Features

- For general laboratory
- Automatic processing
- The control Function of the pressure exhaust speed
- Exact Electric Timer
- Standard Basket
- Four kinds model

Temp. Range	Ambient to 135 °C
Temp. Accuracy	±1 °C
Operating Pressure	1.0 to 1.35kg/cm ²
Temp. Controller	Digital P.I.D. controller
Timer	Electronic Type, 0-60min.
Pressure Gage	Mechanical Type, 0-5kg/cm ²
Air Exhaust	Adjustable Valve
Safety Device	Over Heat protector
Monitor unit	Audible & Visible Device
Basket	Standard Type, 2EA

Model	VS-1221-45	VS-1221-60	VS-1221-80	VS-1221-100
Dimension (In)	Dia.300×630mm	Dia.350×630mm	Dia.400×630mm	Dia.450×630mm
(OUT)	670×470×860(H)mm	670×470×875(H)mm	770×570×910(H)mm	770×570×910(H)mm
Capacity	45 ℓ	60 ℓ	80 ℓ	100 ℓ
Heater	2Kw	3Kw	4Kw	5Kw
Weight	73Kg	81Kg	105Kg	126Kg
Power	220V, 9.5A, 60Hz	220V, 13.6A, 60Hz	220V, 18A, 60Hz	220V, 23A, 60Hz