

Blood Bank Refrigerator

VS-1302MBR2

Temp Control	Micro Processor Controller
Printer	YJ-320 S/T (Thermal type)
Temperature Range	0°C TO 10°C
Temperature Uniformity	±2°C at 4°C
Dimension	(in) 510×525×650(H)mm (out)900×700×740(H)mm
Capacity	174 ℓ
Door	1Door
Compressor	1/5HP, R-134a
Display Type	7-Segment Display
Material	(in) Stainless steel (out)Stainless steel & steel with powder coating
Safety Device	Exclusive Low & High Temp, Protector
Power	220V, 250W, 60Hz
Weight	78Kg



VS-1302SBR

Temp Control	Micro Processor Controller
Printer	YJ-320 S/T (Thermal type)
Temperature Range	0°C TO 10°C
Temperature Uniformity	±2°C at 4°C
Dimension	(in) 620×690×1430(H)mm (out)740×800×2020(H)mm
Capacity	612 ℓ
Door	1Door, Tempered Triple Safety Glass
Compressor	1/3HP, R-134a
Display Type	7-Segment Display
Material	(in) Stainless steel (out)Stainless steel & steel with powder coating
Safety Device	Exclusive Low & High Temp, Protector
Power	220V, 750W, 60Hz
Weight	161Kg



VS-1302LBR

Temp Control	Micro Processor Controller
Printer	YJ-320 S/T (Thermal type)
Temperature Range	0°C TO 10°C
Temperature Uniformity	±2°C at 4°C
Dimension	(in) 1140×690×1430(H)mm (out)1257×800×2020(H)mm
Capacity	1125 ℓ
Door	2Door, Tempered Triple Safety Glass
Compressor	1/3HP, R-134a
Display Type	7-Segment Display
Material	(in) Stainless steel (out)Stainless steel & steel with powder coating
Safety Device	Exclusive Low & High Temp, Protector
Power	220V, 800W, 60Hz
Weight	218Kg