




EASY FITT series (No.10-32UNF, For OD1/16" Tubing)




P/N 9001
EASY FITT φ 10
Carbon-PEEK
Quantity: 10



P/N 9002
EASY FITT φ 12
Carbon-PEEK
Quantity: 10




P/N 9130
EASY FITT TYPE 2
Carbon-PEEK
Quantity: 10




P/N 9004
STOP FITT
Carbon-PEEK
Quantity: 10

HEX FITT series (No.10-32UNF, For OD1/16" Tubing)




P/N 9003
HEX FITT
Carbon-PEEK
For 1/16" tubing
Quantity: 10




P/N 9135
HEX FITT MINI
Carbon-PEEK
For 1/16" tubing
Quantity: 10


DOUBLE LOCK series (No.10-32UNF, For OD1/16" Tubing)



P/N 9301
1/16 DOUBLE-LOCK FITTING S
Stainless Steel 316
Quantity: 10




P/N 9302
1/16 DOUBLE-LOCK FITTING L
Stainless Steel 316
Quantity: 10




P/N 9401
1/16 DOUBLE-LOCK FERRULE
Carbon-PEEK
Quantity: 10


HIGH PRESSURE series (No.10-32UNF, For OD1/16" Tubing)



P/N 9501
1/16 FITTING S
Stainless Steel 316
Quantity: 10




P/N 9502
1/16 FITTING L
Stainless Steel 316
Quantity: 10




P/N 9601
1/16 FERRULE
Stainless Steel 316
Quantity: 10


FLAT SEAL series (No.1/4-28UNF, For OD1/16" Tubing)




P/N 9721
1/16 FLAT SEAL FITTING S
Carbon-PEEK
Quantity: 10



P/N 9703
1/16 FLAT SEAL FERRULE
ETFE
Quantity: 10




P/N 9720
1/16 FLAT SEAL PLUG
PEEK
Quantity: 10




P/N 9720P
1/16 FLAT SEAL PLUG
PP
Quantity: 10

FLAT SEAL series (No.1/4-28UNF, For OD2.00mm Tubing)




P/N 9730
 φ2 FLAT SEAL FITTING
 Carbon-PEEK
 Quantity: 10




P/N 9700
 φ2 FLAT SEAL FERRULE
 ETFE
 Quantity: 10


FLAT SEAL series (No.1/4-28UNF, For OD3.0mm & OD1/8" Tubing)



P/N 9711
 1/8 FLAT SEAL FITTING
 Carbon-PEEK
 Quantity: 10




P/N 9701
 φ3 FLAT SEAL FERRULE
 ETFE
 Quantity: 10




P/N 9702
 1/8 FLAT SEAL FERRULE
 ETFE
 Quantity: 10


FLAT SEAL series (No.1/4-28UNF, For OD1/16" Tubing)



P/N 9748
 1/16 FLAT SEAL FITTING
 Carbon-PEEK
 Quantity: 10




P/N 9742
 1/16 FLAT SEAL FITTING
 POM
 Quantity: 10




P/N 9706
 1/16 M6 FLAT SEAL FERRULE
 ETFE
 Quantity: 10


FLAT SEAL series (No.1/4-28UNF, For OD2.0mm Tubing)



P/N 9747
 φ2 FLAT SEAL FITTING
 Carbon-PEEK
 Quantity: 10




P/N 9741
 φ2 FLAT SEAL FITTING
 POM
 Quantity: 10




P/N 9705
 φ2 M6 FLAT SEAL FERRULE
 ETFE
 Quantity: 10


FLAT SEAL series (No.1/4-28UNF, For OD3.0mm Tubing)



P/N 9746
 φ3 FLAT SEAL FITTING
 Carbon-PEEK
 Quantity: 10




P/N 9740
 φ3 FLAT SEAL FITTING
 POM
 Quantity: 10




P/N 9704
 φ3 M6 FLAT SEAL FERRULE
 ETFE
 Quantity: 10


FLAT SEAL series (M6, For OD1/16" Tubing)



P/N 9753
1/16 M6 FLAT SEAL FITTING
Carbon-PEEK
Quantity: 10




P/N 9734
1/16 M6 FLAT SEAL FITTING
POM
Quantity: 10




P/N 9706
1/16 M6 FLAT SEAL FERRULE
ETFE
Quantity: 10


FLAT SEAL series (M6, For OD2.0mm Tubing)




P/N 9752
2 M6 FLAT SEAL FITTING
Carbon-PEEK
Quantity: 10



P/N 9732
2 M6 FLAT SEAL FITTING
POM
Quantity: 10




P/N 9762
2 M6 FLAT SEAL FITTING
PP
Quantity: 10




P/N 9705
2 M6 FLAT SEAL FERRULE
ETFE
Quantity: 10


FLAT SEAL series (M6, For OD3.0mm Tubing)




P/N 9751
3 M6 FLAT SEAL FITTING
Carbon-PEEK
Quantity: 10



P/N 9731
3 M6 FLAT SEAL FITTING
POM
Quantity: 10




P/N 9761
3 M6 FLAT SEAL FITTING
PP
Quantity: 10




P/N 9704
3 M6 FLAT SEAL FERRULE
ETFE
Quantity: 10

FLAT SEAL series (M6, For OD1/8" Tubing)




P/N 9760
1/8 M6 FLAT SEAL FITTING
PP
Quantity: 10




P/N 9707
1/8 M6 FLAT SEAL FERRULE
ETFE
Quantity: 10


DOUBLE LOCK series (M8, For OD3.0mm & OD1/8" Tubing)




P/N 9304
1/8 DOUBLE LOCK FITTING
Stainless Steel 316
Quantity: 10




P/N 9305
1/8 DOUBLE LOCK FITTING
PP
Quantity: 10



P/N 9402
1/8 DOUBLE LOCK FERRULE
ETFE
Quantity: 10




P/N 9405
3 DOUBLE LOCK FERRULE
ETFE
Quantity: 10




P/N 9407
1/8 DOUBLE LOCK FERRULE
PEEK
Quantity: 10


UNION




P/N 9005
UNION
Carbon-PEEK
No.10-32UNF
Quantity: 1




P/N 9511
UNION
Stainless Steel 316
No.10-32UNF
Quantity: 1



P/N 9736
M6 FLAT SEAL UNION
POM
M6
Quantity: 1




P/N 9101
1/4-28 FLAT SEAL UNION
Carbon-PEEK
M6
Quantity: 1




P/N 9111
UNION
Carbon-PEEK
M8
Quantity: 1


REDUCING UNION




P/N 9110
REDUCING UNION
Carbon-PEEK
No.10-32UNF, No.1/4-28UNF
Quantity: 1



P/N 9008
REDUCING UNION
Carbon-PEEK
M8, No.10-32UNF
Quantity: 1




P/N 9737
M6-10-32 FLAT SEAL UNION
POM
M6, No.10-32UNF
Quantity: 1




P/N 9112
SYRINGE UNION
Carbon-PEEK
No.10-32UNF, Luer tip syringe connector
Quantity: 1


TEE




P/N 9006
TEE
Carbon-PEEK
No.10-32UNF
Quantity: 1



P/N 9131
1/4-28 FLAT SEAL TEE
Carbon-PEEK
No.1/4-28UNF
Quantity: 1




P/N 9121
TEE
Carbon-PEEK
M8
Quantity: 1




P/N 9738
M6 TEE
POM
M6
Quantity: 1

REDUCING TEE




P/N 9108
REDUCING TEE
Carbon-PEEK
No.1/4-28UNF, No.10-32UNF
Quantity: 1




P/N 9107
REDUCING TEE
Carbon-PEEK
M8, No.10-32UNF
Quantity: 1

CROSS



P/N 9007
CROSS
Carbon-PEEK
No.10-32UNF
Quantity: 1



P/N 9109
CROSS
Carbon-PEEK
No.1/4-28UNF
Quantity: 1

PEEK Tubing



P/N 5999
Natural PEEK
Red Striped
OD 1/16"
ID 0.13mm
Length: 10m



P/N 6000
Natural PEEK
Yellow Striped
OD 1/16"
ID 0.17mm
Length: 10m



P/N 6001
Natural PEEK
Blue Striped
OD 1/16"
ID 0.25mm
Length: 10m




P/N 6002
Natural PEEK
Orange Striped
OD 1/16"
ID 0.50mm
Length: 10m




P/N 6003
Natural PEEK
Green Striped
OD 1/16"
ID 0.75mm
Length: 10m



P/N 60031
Natural PEEK
Gray Striped
OD 1/16"
ID 1.0mm
Length: 10m




P/N 6004
Natural PEEK
No Stripe
OD 1/8"
ID 1/16"
Length: 5m



P/N 60042
Natural PEEK
No Stripe
OD 1/8"
ID 2.0mm
Length: 5m

PTFE Tubing



P/N 7005
PTFE
OD 3.0mm
ID 2.0mm
Length: 10m



P/N 7211
PTFE
OD 1/16"
ID 0.25mm
Length: 10m



P/N 7210
PTFE
OD 1/16"
ID 0.5mm
Length: 10m



P/N 7212
PTFE
OD 1/16"
ID 0.8mm
Length: 10m

ETFE Tubing



P/N 7009
ETFE
OD 1/16"
ID 0.2mm
Length: 10m



P/N 7010
ETFE
OD 1/16"
ID 0.3mm
Length: 10m




P/N 7011
ETFE
OD 1/16"
ID 0.5mm
Length: 10m



P/N 7012
ETFE
OD 1/16"
ID 0.8mm
Length: 10m




P/N 7013
ETFE
OD 1/16"
ID 1.0mm
Length: 10m




P/N 7016
ETFE
OD 2.0mm
ID 1.0mm
Length: 10m


PFA Tubing



P/N 7103
PFA
OD 2.0mm
ID 1.0mm
Length: 10m




P/N 7104
PFA
OD 3.0mm
ID 2.00mm
Length: 10m




P/N 7105
PFA
OD 1/8"
ID 1/16"
Length: 10m

Solvent Filter



P/N 8800
Suction Filter
PP
Pore size 5 μ m
For ID2.0mm tubing
Quantity: 1



P/N 8802
Suction Filter
PTFE, Stainless Steel
Pore size 10 μ m
For ID2.0mm tubing
Quantity: 1