

Swing Rotor, VR-50S



Available unit : VS-5500CFi, VS-5500N, VS5000N


Max. RPM	Max.RCF	Rotor Capacity	Radius
5500RPM	5580×g	480ml, 15ml×32	167mm

Bucket



Model	VB-32	VB-08	VB-16	VB-04
Tube	15ml	50ml	15ml	50ml
Places	8	2	4	1
Tube Dim.φ×L mm	16×114	29×104	17×120	30×115
Max. RPM / VS-5500N	4,500 RPM	5,000 RPM	4,500 RPM	5,500 RPM
VS-5000N	4,000 RPM	4,500 RPM	4,000 RPM	5,000 RPM

Rack Carrier



Model	VB-40	Model	RK-32	RK-48	RK-56	RK-80
Max. RPM	3500RPM	Tube	15ml	10ml	5ml	3ml
		Places	8	12	14	20
		Total(Tubes)	32	48	56	80

Micro Plate Rotor



Model	VR-02M
Max. RPM	2,500RPM
Max. RCF	2,000×g
Capacity	96well Micro plate×2EA