


HPLC New

Belmont s/ Yverdon Thursday, 10 February 2011




The easiest way to work with your HPLC, by using Clarity software from Data Apex.


Configuration in 3 steps only :



Choose your instruments included in your system
Insert them into the right side window and valid the configuration



Define the different instrument parameter
Define the method, flow rate, gradient shape, detector wavelength, sensitivity, Autosampler protocol,
Create your working list
Define your data treatment and store it in your application file



Start your application, start the HPLC Pump / GC, and run your samples until you get the results on your display.
Check your results on the data treatment display, optimize the integration, print out results It's done...!

If you need other detector, like DAD, or use other brand name instruments together, use the same 3 step configuration !

Price start from 4200 CHF ex VAT / MWST 8%

For future information's, demo CD or quote contact us at :

TecDev
Le Villaret
CH -1432 Belmont s/Yverdon
Schweiz
Email: tecdev@hispeed.ch Tel: +41 (0) 24 435 22.02

Tec Dev Le Villaret 1432 Belmont s/ Yverdon Switzerland Phone : ++41 (0) 24 435 2202 Fax : ++41 (0) 24 435 2203 E mail : tecdev@swissonline.ch

Banking details UBS AG Yverdon Account No : G 3.323.123.1
Swift code : UBS Yverdon les Bains : UBSWCHZH14A
IBAN code : CH 45 0029 7297 G 332 3123
VAT / TVA / MWST No : 553 253

Managing director :
Jean Pierre Bloch