

Small sized High Speed Refrigerated Centrifuge

Model: VS-15CFi



Feature


Max. 17,000 RPM, Max. 26,171 xg
Inverter Motor Drive
Full Microprocessor control
Safety Apparatus
Easy Loading

Max. RPM	17,000 RPM
Max. RCF	26,171 × g (Angle rotor)
Max. Capacity	600ml, 100ml × 6Tubes
Display	7-Segment Numeric Display
Operating mode	RPM, RCF
Program memory	20 Memory
Drive System	Inverter Motor drive
Acceleration	Slow/ Fast
Deceleration	Slow/ Fast/ Normal
Timer	99Hour. 59Min & Hold for long run
Pre-Condition	Pre-cooling
Temperature Range	-20 ~ 60°C
Safety system	Over speed detection, Unbalance control, Temperature Limit control, Door detector
Dimension	772 × 565 × 848(H)mm
Weight	135kg
Power	220V, 1.6Kw, 50/60Hz
Standard Items	Body, Fixed Angle Rotor V506A, Instruction Manual

Rotor



V1.524A



V1512A



V508A