

Ice Maker (Snow Type)



VS-625NS



VS-625N



VS-625ND



Optional Accessories
Water Filter

Model	VS-625NS	VS-625N	VS-625ND
Production/day	50kg	200kg	400kg
Dimension	(Out)	360×475×615(H)mm	1030×570×1250(H)mm
	(In)	300×410×200(H)mm	425×450×450(H)mm
Compressor	3/8HP, R404a	1HP, R22	2HP
Power Consumption	770W	1.3Kw	1.8Kw
Power	220V, 3.5A, 50/60Hz	220V, 6A, 50/60Hz	220V, 50/60Hz
Weight	33Kg	87Kg	172Kg

Ice Maker (Flake Type)

Model: VS-625



Production	250 Liter/Day
Ice Type	Flake Type
Storage Capacity	159 Liter
Compressor	1/3 HP
Door	Front sliding door
Material (IN)	Stainless steel
	(OUT)
Safety Device	Power breaker, Water level monitor
	Inventory level monitor
Dimension (IN)	500×490×650(H)mm
	(OUT)
Weight	150Kg
Power	220V, 6.7A, 1.3Kw, 50/ 60Hz
Standard Items	Tubing Kit, Scoop