

Vortex Mixer

Model KMC-1300V



Vortex Mixer, Platform Head



Vortex Mixer, Mixing Cup Head

Vortex Mixer

KMC-1300V

Speed Range / Vitesse / Geschwindigkeit

0 - 3000 RPM

Operation Mode / Mode de fonction / Funktion

Touch On & Constant On

Material / Matériel du boîtier

Aluminum Die-casting & ABS

Head / Type de tête / Support

Rubber Mixing Cup

Dimension / Dimensions / Abmessungen

130 x 155 x 155(H) mm, 3.1kg

Power / Alimentation / Speise

AC230V, 0.14A, 60Hz

233.-

Optional Accessory

V130001	Mixing Cup
V130002	Platform Head
V130003	Micro plate Kit
V130004	Multi sample plate Kit



Micro plate Kit



Multi sample plate Kit

TecDev Le Villaret Belmont s/Yverdon Switzerland / Phone +41 (0) 24 435 22 02 Fax: +41 (0) 24 435 22 03 Email tecdev@hispeed.ch

Prix hors TVA 7,6% et transport

Price without MWST 7,6% and shipping cost

Preis ohne MWST 7,6% und transport

Valid until 28 February 2007, Validité de l'action 28 février 2007, Preise Gültig bis 28 Februar 2007