

# Spin Boy

Model VS-100B-M  
VS-100B-N



VS-100B-M

VS-100B-N



- ✚ Mini style
- ✚ Model VS-100 B-M also for vortexer
- ✚ Up to 8 PCR tube strip
- ✚ Quick spin down
- ✚ 8 holes Rotor
- ✚ Funny Neon light
- ✚ Funny price !

## Multi Mini Centrifuge VS-100B-M

Speed	10,000 RPM / 4000RPM
G-force	5,590 x g
Function	1. Centrifuge: 1.5ml x 8 Tubes 2. Centrifuge: 0.2ml 8 PCR Tube Strips x 2 3. Vortex Mixer
Drive	Direct drive system
Dimension	153 x 162 x 122(H) mm
Lighting	Neon
Power	220V, 60Hz

**Configuration:** main body, 1,5 ml x 8 ml tube rotor Vortex head, PCR tube rotor

## Mini Centrifuge VS-100B-N

Speed	10,000 RPM
G-force	5,590 x g
Capacity	1.5ml x 8 Tubes
Drive	Direct drive system
Dimension	153 x 162 x 122(H)mm
Lighting	Neon
Power	220V, 60Hz

**Configuration:** main body, 1,5 ml x 8 ml tube rotor