


Fridge XXL™

Sample cooling rack



- Small profile cooled rack for large sample volume
- Fit in all standard Robotic system or as stand alone module
- Compatible with any standard samples vials or tubes
- Cooling range :
4 °C to RT+5 °C standard version
4 °C to 60 °C extended version
- Can be used with detachable rack
- 175 x 240 mm rack for 154 samples tube of 11 x 70 mm
- 240 x 240 mm rack for 100 samples tube of 16 x 10 mm
- Different rack available
- Use external temperature controller with PID function
- Fast temperature stabilisation and high temperature precision of up to 0,1 °C
- Base plate dimension : 545 x 340 mm.



<u>Reference</u>	<u>Description</u>	<u>Price (US \$)</u>
U Fridge 711-XXL	U Fridge XL cooling rack 320 W version	2470.-
U Fridge 711-XXL ext	U Fridge XL cooling / heating rack 320 W extended version	2970.-
KU-EU-XXL	PID Controller base size version complete 115 V/ 230 V	2270.-
KU-EU-XXL ext	PID Controller extended size version complete 115 V/ 230 V	2570.-
REM	Remote option	250.-
R154 x 11	Detachable rack 154x 11 mm diameter holes	440.-
R100 x 16	Detachable rack 100 x 16 mm diameter holes	420.-
R6 x MTP	Detachable rack 6 x 96 / 384 / 1536 vials microtiterplate	310.-
R25x24	Detachable rack 25 x 24 mm diameter scintillation vials	380.-
R154x15-Epp	Detachable rack 154 x 10 mm Eppendorf 1,5 ml	440.-
R300x6	Detachable rack 300 x 6 mm Eppendorf 0,6 ml	610.-
R-CUST	Customised rack	ask
Logo	Your Logo on front plate of controller	50.-

Delivery: 5-6 Week
Warranty: 1 year
Price FOB: Belmont s/ Yverdon

Copyright TecDev 6 2002

TecDev
Le Villaret
1432 Belmont s/Yverdon
Switzerland

Tel : ++41 (0) 24 435 2202
Fax : ++41 (0) 24 435 2203
E mail : tecdev@swissonline.ch