


Fridge™



- Small cooled rack for reagent or samples cooling
- Fit in all standard laboratory sampler or laboratory robot
- Small foot print (w x l x h) : 100 x 75 x 100 mm
- Cooling range :
 - 7 °C to RT-5 °C standard version
 - 7 °C to 40 °C extended version
 - (RT = 20 degree)
- Use detachable rack with different holes diameter
- Fit with any standard vessels
- Cooling volume of 88 x 53 x 30 mm
- Use external temperature controller with PID function
- Allow fast temperature stabilisation and high temperature precision of up to 0,1 °C.



<u>Reference</u>	<u>Description</u>	<u>Price (US\$)</u>
U-fridge	U Fridge rack standard version	520.-
U-fridge-XL	U Fridge rack extended version	570.-
KU-EU	PID Controller standard version complete 115 V/ 230 V	950.-
KU-EU XL	PID Controller extended version complete 115 V/ 230 V	1150.-
INS-2x26	Detachable rack 2 x 26 mm diameter holes	170.-
INS-8x12	Detachable rack 8 x 12 mm diameter holes	220.-
INS-5x16	Detachable rack 5 x 16 mm diameter holes	220.-
INS-30x6	Detachable rack 30 x 6 mm diameter holes	310.-
INS-12x8	Detachable rack 12 x 8 mm Eppendorf 1,5 ml	250.-
INS-24x6	Detachable rack 30 x 6 mm Eppendorf 0,6 ml	310.-
INS-CUST	Customised rack	320.-
Logo	Your Logo on front plate of controller	50.-

Delivery: 4-5 Weeks

Warranty 1 year

Price FOB: Belmont s/ Yverdon

TecDev

Le Villaret

CH 1432 Belmont s/ Yverdon

Email : tecdev@swissonline.ch

Tel : ++41 (0) 24 435 2202

Fax : ++41 (0) 24 435 2203